

# Stop reagent • Instructions for use • IVD CE

REF: RLR011, RLA703

## Intended purpose & product description

The device is a consumable, which is used for performing the basophil histamine release assay (BHRA).

The stop reagent is included the assay test completion, and serves as a stabilizer in the coupling procedure by lowering the pH.

For in vitro diagnostic use.

## Intended user

The stop reagent is for professional (health care) use only.

## Composition

Perchloric acid 0.59 % w/v (0.059 mol/L)

## Variants

Product name	Article no	UDI code	Volume	Pack size
Stop reagent, 1 L	RLA703	05744003280254	1 L	1
Stop reagent, 100 mL	RLR011	05744003280247	100 mL	1 (as part of the RLA700 Reagent Kit)

## Reagent preparation

Please refer to the Instructions for use for the OPA disc (Art. no. RLR013, RLA708 and RLA708X) for preparation of the coupling reagent and description of the coupling procedure.

## Test principle

Refer to separately provided BHRA / HR-Test method-specific work instructions.

## Storage and handling conditions

Store at +2 - +8 °C until labelled expiry.

The reagent should reach room temperature before use.

## Warnings

The stop reagent is a non-hazardous reagent

## Waste disposal

Non-hazardous waste










# Stop reagent • Instructions for use •

REF: RLR011, RLA703

## Manufacturer

RefLab ApS  
Copenhagen Bioscience Park 'COBIS'  
Ole Maaloes Vej 3  
DK 2200 Copenhagen N  
Denmark

## Symbols

	CE mark
	For in vitro diagnostic use
	Catalog number
	Batch code
	GTIN code
	Expiry date
	Storage temperature
	Consult instructions for use
	Manufacturer