

## Requisition of HR-Urticaria Test



Member of R-Biopharm Group

RefLab • P.O. BOX 590 • Ole Maaløes Vej 3 • DK- 2200 Copenhagen N  
Phone: +45 7070 2345 • Fax +45 7070 2367 • [reflab@reflab.dk](mailto:reflab@reflab.dk) • [www.reflab.dk](http://www.reflab.dk)

PATIENT	CLINIC
<b>Name</b> _____	_____
<b>Address</b> _____ _____	_____
_____	_____
<b>Date of birth</b> _____	_____
<b>Date of blood sampling</b> _____	_____
_____	_____
<b>Male</b> <input type="checkbox"/> <b>Female</b> <input type="checkbox"/>	

Detection of autoantibodies against the IgE and/ or IgE receptor

HR-Urticaria Test

### Practical information – Please note

Serum (**1-2 ml per patient placed in a transport container**) is posted by ordinary mail at ambient temperature to RefLab.

**Method:** Patient serum and donor basophils are incubated. Histamine is measured in stimulated and un-stimulated cells and the results are expressed as histamine release in percentage of the cellular histamine content.

A histamine release > 16,5 % is a positive test result.

The test result will be returned to the clinic within 10 days after receipt to the laboratory.