

Requisition of HR-CSU Test



RefLab • Ole Maaløes Vej 3 • DK- 2200 Copenhagen N • Denmark
Phone: +45 7070 2345 • Fax +45 7070 2367 • info@reflab.dk • www.reflab.dk

PATIENT	CLINIC
Name
Date of birth	
Date of blood sampling	Email:
Male <input type="checkbox"/> Female <input type="checkbox"/>	

Detection of autoantibodies against the IgE and/ or IgE receptor

HR-CSU Test

Practical information – Please note

Serum (**1-2 ml per patient placed in a transport container**) is send at ambient temperature to RefLab.

Method: Patient serum and donor basophils are incubated. Histamine is measured in stimulated and un-stimulated cells and the results are expressed as histamine release in percentage of the cellular histamine content.

A histamine release > 16,5 % is a positive test result.

The test result will be returned to the clinic within 10 days after receipt to the laboratory.

Please fill in contact details on where to send the test result.