

Requisition of s-BHRA/HR-Urticaria Test



Member of R-Biopharm Group

RefLab ApS • Ole Maaløes Vej 3 • 2200 Kbh. N • Tlf: 70 70 23 45 • Fax: 70 70 23 67 • info@reflab.dk • www.reflab.dk

PATIENT	CLINIC
Name:
Address:
Date of birth:	BILLING INFORMATION
Sample collection date:

Detection of auto-antibodies against IgE and/ or the IgE-receptor
s-BHRA/HR-Urticaria Test

Practical information – Please note:

Serum (1-2 mL per patient placed in a transport container) is posted by ordinary mail at ambient temperature to RefLab ApS.

Method: Patient serum and donor basophils are incubated. Histamine is measured from stimulated and un-stimulated cells and the results are expressed as histamine release in percentage compared to the total cellular histamine content.

A histamine release >16.5 % is a positive test result.

The test result will be returned to the clinic within 10 days after receipt to the laboratory.