

Requisition of s-BHRA/HR-Urticaria Test



Member of R-Biopharm Group

RefLab ApS • Kanalholmen 1, 2.th. • 2650 Hvidovre • Denmark • Tlf: 70 70 23 45 • info@reflab.dk • www.reflab.dk

PATIENT	CLINIC
Name:
Address:
.....
.....	BILLING INFORMATION
Date of birth:
Sample number:
Sample collection date:
PO-number/ Cost Center:

Detection of auto-antibodies against IgE and/ or the IgE-receptor

s-BHRA/HR-Urticaria Test

Practical information – Please note:

Serum (1-2 mL per patient placed in a transport container) is posted by ordinary mail at ambient temperature to RefLab ApS.

Method: Patient serum and donor basophils are incubated. Histamine is measured from stimulated and un-stimulated cells and the results are expressed as histamine release in percentage compared to the total cellular histamine content.

A histamine release >16.5 % is a positive test result.

The test result will be returned to the clinic within 7-21 business days after receipt to the laboratory.